

fluorescent in the wrong parameter, it can be found before the experimental tube is run). Do this process for each control tube. It is critical that PMT voltages are not changed once you begin adjusting compensation. If it becomes necessary, all the controls must be run again to ensure compensation is correct. In addition, further tubes are required to properly set the compensation for the CD45/FITC and CD14/RPE-Cy5 parameter. In this case, samples stained with CD3/FITC and CD3/RPE-Cy5 can be used as the compensation controls. Hint: Once you have set the compensation for all parameters, save the protocol you are using in Summit Software. This way, next time you repeat your experiment (using the same antibodies) the compensation will be just about perfect and may only need a little touch up. Once the compensation is complete you are ready to run the experimental sample which is shown in Figures 2-8. Populations were labeled using the annotation feature. Color gating was applied to highlight distinct populations and/or make the data more presentable. The workspace can be exported as graphics as a whole or individual histograms/plots and statistics to Microsoft® Excel.

Discussion

As shown, multicolor immunophenotyping is becoming much easier with the new instrumentation available. By using more markers in each panel, fewer tubes are required using fewer sample and less time. The CyAn ADP LX 9 Color offers an excellent platform for multicolor analysis with the small

bench-top footprint, fixed alignment, spatially separated lasers and easily interchangeable filters. The software allows the user to save the PMT settings, compensation, and gating which makes repetition of an experiment easy. Compiling data into a report format is simple with cut and paste tools readily available.

TECHNICAL TIPS

- > During cell preparation single controls are used in this protocol to demonstrate instrument set up for a given experiment. Although single controls are always a good idea to troubleshoot staining problems, once the experiment has been set up there is very little variation run to run. The blood used in the cell preparation should not be older than 24 hours. When using the Uti-Lyse, Reagent A is the fixative. Therefore, the 10 minutes incubation time for reagent A must be adhered to so as to avoid fixing red cells as well as the white cells. Reagent B is the lysing agent. As the wash steps occur after the 10 minutes lyse, very few red cells should be in the pellet after centrifugation. Cells can be stored at 4 °C / 39.2 °F for up to 24 hours before analysis. After 24 hours, the scatter properties will begin to deteriorate.
- > Before running samples verify that the instrument has had the appropriate time to warm up and the alignment specifications have been met using a standard bead particle.
- > During acquisition single controls should be run to properly set voltages and compensation. Once set, the protocol can be saved which will contain all the settings and can be reused for this same experiment in the future.
- > During analysis color gating can be done in Summit Software to highlight populations and show marker overlap.

PRODUCT	CODE
CyAn ADP LX 9 Color	CY201
Monoclonal Mouse Anti-Human CD3/FITC	F0818
Monoclonal Mouse Anti-Human CD3/RPE-Cy5 ..	C7067
Monoclonal Mouse Anti-Human CD45/FITC	F0861
Monoclonal Mouse Anti-Human CD56/RPE	R7251
Monoclonal Mouse Anti-Human CD19/APC	C7224
Uti-Lyse, Erythrocyte Lysing Reagent	S3350
Hank's Buffered Salt Solution	S1754

For research use only – not to be used in diagnostic procedures. Other vendor products used in this application: eBioSciences and Caltag/Invitrogen.

The protocols in this application note might deviate from the normal recommended protocol/specification guidelines which are included with the Dako product or any other non-Dako product.

Graphs

Figures 1 to 8

Gating strategy:

Fig. 1: Create "Lymphocyte gate for compensation".

Fig. 2: CD3/SSC to create "T Lymphocyte gate"

Fig. 3: CD4/CD8 plot using the combination of "Final Lymphocyte gate" and the "T Lymphocyte gate" to identify CD4+ and CD8+ T Lymphocytes.

Fig. 4: CD45/SSC plot showing initial lymphocytes gate ("CD45 vs. SSC Lymphocyte gate").

Fig. 5: CD45/CD14 plot using CD45 vs. SSC Lymphocyte gate to create the "Final Lymphocyte gate", which excludes monocyte contamination.

Fig. 6: CD56/CD8 plot using "Final Lymphocyte gate" to identify the presence of CD8+/CD56+ cells

Fig. 7: CD3/CD56 plot using "Final Lymphocyte gate" to identify NK Lymphocytes (CD56+/CD3-) and NK like Lymphocytes (CD3+/CD56+).

Fig. 8: CD3/CD19 plot using "Final Lymphocyte gate" to identify B lymphocytes.

Figure 1




Figure 2




Figure 3




Figure 4




Figure 5




Figure 6




Figure 7




Figure 8

