




*Better Tools. Good Science.*



**The Summit CICERO Upgrade improves the performance of flow cytometers by replacing the outmoded computer hardware, software and electronics systems with Cytomation's Summit Software System, Summit Workstation and CICERO Electronics Interface.**

## **SUMMIT SOFTWARE SYSTEM**

Summit's software compensation allows you to display raw or compensated values, and interactively adjust constants directly on the histogram or through the Compensation Matrix. Compensation values can be displayed on the axis, so the constants are updated in real-time and are displayed on graphs when printed. Compensation parameters can be defined for an unlimited number of colors. Compensated and raw data can be displayed side-by-side to facilitate the creation of "what-if" scenarios. Constants can be re-used across samples, or separate compensation matrices can be created for each of an unlimited number of samples, all on the screen at the same time. Compensation is available in a 3 x 3 matrix (including FL1:FL3), ratios and time.

Summit offers flexibility throughout the program so that data can be gathered and manipulated in ways that best suit the application and user. For example, a logical expression can be defined using any number of regions in "AND" and "NOT" combinations, including non-rectangular regions with an unlimited number of vertices. Up to 32 analysis and gating regions can be specified for each histogram.

Summit's Data Navigator tracks samples in experiments, projects, and research programs. Tests and controls with different parameters can be grouped, and samples can be dragged and dropped between experiments to help organize information. The desktop can be organized by editing tables, changing histogram size preferences and moving around toolbar buttons. Each user's desktop organization preferences can be saved. All organizational tools available in Windows® NT are available in Summit.

Summit graphics, statistics, raw data, compensation matrices, and other useful information are easily copied to the clipboard and pasted into other applications. Summit supports color graphics output on any printer available in Windows NT. The ability to customize fonts, colors, and sizes provides desktop publishing quality output.

## **SUMMIT WORKSTATION**

This high-end computer workstation provides a powerful platform for the Summit CICERO Upgrade, as well as off-line analysis, manuscript preparation and desktop publishing. The standard workstation is a very powerful, network-ready computer, and upgrades and peripherals are available to customize each workstation to meet the needs of the researcher.

## CICERO ELECTRONICS INTERFACE

The CICERO Electronics Unit interfaces with most common flow cytometers and improves their performance. It's high-speed data acquisition, in excess of 20,000 events per second, reduces data loss when studying rare events.

The CICERO Electronics Interface has 16 bit ADCs (analog-to-digital converters), which provide high-resolution data capture.

## SUMMIT CICERO TECHNICAL SPECIFICATIONS

### Summit CICERO Upgrade includes:

Summit Software CD-ROM with program files  
Summit Workstation (see specifications below)  
CICERO Electronics Interface with interface cables  
Summit Getting Started Guide  
Summit Reference Guide

### Summit Software Platform:

Windows® NT for acquisition, Windows NT, 98, or 95 for offline data analysis

### Optical Specifications:

Determined by flow cytometer optical system

### Fluidics Specifications:

Determined by flow cytometer fluidics system

## SUMMIT WORKSTATION SPECIFICATIONS

### Platform:

Windows NT operating system

### Motherboard:

Dual Intel 400 MHz Pentium® II with 128 MB RAM expandable to 1 GB

### Storage Space:

Two 8.4 GB hard disks, CD-ROM drive, optional DVD, CD-RW

### Monitor:

Single or dual-monitor 16 MB video card, one or two 21-inch or larger traditional tube-type or LCD flat screen monitors

### Miscellaneous:

Ethernet LAN card, 3.5" disk drive, sound card and speakers

Summit CICERO Upgrades are available for analysis work on the following systems. Ask your sales representative for information about Cytomation's Summit MoFlo Electronics Upgrade for instruments performing sorts.

- FACScan
- EPICS ELITE
- EPICS C
- FACStar+
- EPICS 700 Series
- Ortho 50
- FacStar
- EPICS 500 Series
- Partec
- FACS IV
- EPICS Profile

Cytomation, Inc. is a privately held bio-technical instrumentation corporation specializing in high-performance, high-speed flow cytometer analyzers, sorters and upgrades. Our mission is to design, produce, and service the finest flow cytometers and cell sorters in the world — unparalleled in performance, accuracy, versatility, reliability and speed. MoFlo, our premier flow cytometer, is a modular system that is easily upgraded whenever requirements change or new modules become available.

### North America

Cytomation, Inc.  
4850 Innovation Drive  
Fort Collins, CO 80525  
USA  
Telephone: 970.226.2200  
Facsimile: 970.226.0107  
E-mail: [info-northamerica@cytotation.com](mailto:info-northamerica@cytotation.com)

### Europe

Cytomation Bioinstruments GmbH  
Engesserstrasse-4  
79108 Freiburg im Breisgau  
Germany  
Telephone: 49.761.559.7370  
Facsimile: 49.761.559.7375  
E-mail: [info-europe@cytotation.com](mailto:info-europe@cytotation.com)

### Australasia

Cytomation PTY LTD  
60-66 Hanover Street  
Fitzroy, Vic 3065  
Australia  
Telephone: 613.9417.4577  
Facsimile: 613.9417.2987  
E-mail: [info-australasia@cytotation.com](mailto:info-australasia@cytotation.com)